

# Working in the Blue Sky Plan Planning Software


## A Quick Guide for Thommen Medical Users

The guiding sleeve for VECTOdrill™ pilot drill from Thommen Medical can now be selected for planning with Blue Sky Plan:


**1**  
In the first planning step, one of three pilot drill length options needs to be selected.

### Available drill lengths


Thommen Pilot Drill **extra short** –  
Thommen Medical Sleeve



Thommen Pilot Drill **short** –  
Thommen Medical Sleeve




Thommen Pilot Drill **long** –  
Thommen Medical Sleeve



**2**  
Based on the selected drill length, certain implant lengths are available for the detailed case planning in the second step.

### Implant length (mm)

6.5    8.0    9.5    11.0    12.5    14.0    17.0



Visit our website [thommenmedical.com](http://thommenmedical.com) for more information on guided surgery solutions or contact us via [tmdigital@thommenmedical.com](mailto:tmdigital@thommenmedical.com)